## PULSTRODE PLUS

REF PT46-4

## EASYTRODE PLUS

REF ET46-4

# MULTIPOLAR ELECTROCATETER FOR RADIOFREQUENCY NEUROMODULATION TREATMENT AND EPIDURAL ADHESIONS LYSIS

Please read carefully the instructions for a proper use of the device.

### DESCRIPTION

These devices are small diameter leads specifically designed for spinal epidural space applications. Once introduced into space, they can be used for pharmacological, mechanical (lysis) and electrical treatments. For electrical applications of this device, it is necessary the connection to an external equipment that provides a Pulsed Radiofrequency signal.

### INDICATIONS

Epidural invasive diagnosis and treatments for acute and chronic pain through centralized epidural treatments, such as:

- Mechanical by making lysis of adhesions and/or stenosis of the vertebral epidural canal, due to degenerative pathophysiological processes or cicatricial outcomes from other surgery interventions;
- Pharmacological by administering injections directly to the epidural space, for all pathologies curable by drug administration directly in epidural space target;
- Electrical by performing pulsed RF @ 42° Celsius degree treatments with the use of a special external equipment that provides compatible signal;
- Low frequency electrical, for stimulation and neuromodulation of epidural nerve structures for the study of the effectiveness of a totally implantable SCS system.

### CONTRAINDICATIONS

These devices are contraindicated when, in the opinion of the physician, they may be contrary to the patient's health. In particular in patients with coagulation disorders, in the treatment of anticoagulants; patients with high cardiovascular risk; patients with ongoing infections; patients who are unable to hold their position; patients with psychiatric disease, pregnant patients and patients under 18.

**CAUTION: Total anesthesia is contraindicated for this type of procedure.**

### INSTRUCTIONS FOR PROPER USE AND WARNINGS

- Prepare the sterile field on the patient body by only discovering the access site.
- Inject a local anesthetic in the access area for next insertion of the steel introducer.
- For both devices the access is via Sacral Hiatus, or intervertebral paramedian to any lumbar, thoracic or cervical region. So explained, make the extradural access with the metal introducer contained in the kit, taking care to don't damage the underlying nervous structures. Use for this operation a low resistance syringe. For PULSTRODE PLUS and EASYTRODE PLUS kit, the low resistance syringe is already included with the others supplied accessories. Therefore connect its luer lock port with introducer. Then proceed under fluoroscopy to position the introducer needle with the right inclination through the interested structures.
- The PULSTRODE PLUS and EASYTRODE PLUS kit, have also included in the Torquer that facilitates the catheter grip by the user. For these PLUS models, insert the tip catheter through the Torquer from right direction (Fig. 1) and then place this in the proximal position of the catheter (Fig. 2).


Fig. 1: Torquer insertion

Fig. 2: Torquer proximal placement

- Under fluoroscopy, insert the catheter through the introducer and reach the extradural epidural space. Place the lead tip with the electrodes in the area affected by the treatment. To facilitate the placing procedure, you can slightly bend the catheter tip near the electrodes by using only your fingers (fig. 3) and perform this manipulation only with the stylet inserted.


Fig. 3: tip bending


Fig. 4: Torquer use

- During the placing of the catheters, handle them by acting on the Torquer and grasping it from the knurled side (fig.4). The Torquer can be positioned anywhere on the outside portion of the catheter, depending on the operator's convenience. To move the Torquer from one position to another, press the white button and while holding down reposition the Torquer and release the button.
- To perform mechanical lysis into the epidural space, in the regions of greatest interest, repeat the movements of the tip in the same area in order to release the epidural space from any fibrotic or stenosis bonds and adhesion which prevent the normal crossing of the catheter in this zone.
- For Pulsed RF @ 42° Celsius degree treatments, now connect the catheter to the special cable connected to the external generator (fig. 6). In order to avoid any contamination of the sterile field, the connecting cable must be inserted beforehand inside a wire-drape. The PULSTRODE PLUS and EASYTRODE PLUS kit already include the plastic protection for covering the cable (fig. 5). Taking care when connecting the two connectors, to close the distal end of the bag with the supplied adhesive tape.

Fig. 5: cable cover use


Fig. 6: catheter - cable connection

- Proceed now with the survey, by means of low frequency stimulation (sensory and/or motor), for the detection of areas to be treated. At this stage, an active and lucid collaboration of the patient is indispensable, which will have to relate the perceptions experienced during stimulation.
- Optimize the position of the catheter in the epidural space, until reaching the overlapping between areas of stimulus perception with the typically painful areas of the patient. For a proper procedure, the threshold voltage value above which the patient begins to perceive sensory stimulation must remain around 0.3 V and not exceed the maximum value of 0.8-0.9 V. For higher values The treatment may be ineffective.
- Identified therefore the optimal area in terms of threshold value and perception area, now perform Pulsed RF @ 42° Celsius degree treatment for a maximum duration of 5 min and voltage as much as max tolerable for the patient, but not exceeding 75 V. It is possible repeat more treatments in different areas during the same operating session. Never run procedures at temperatures above 42° Celsius degree. The sensations experienced by the patient during the treatment may be a pulsation in the target areas of the treatment, tingling, muscle contractions or nothing. Lower the voltage of the treatment or suspend it if any sensations experienced by the patient did not prove to be tolerable.
- Once the treatment has been performed, if it is considered effective by the treating physician, it is possible to remove the stylet from the catheter and infuse through the Luer-Lok port, drugs to complete the procedure. To reinsert the stylet inside the catheter, pull the lead out until get a rectilinear shape and only now reinsert the stylet.
- Both catheters give off a bipolar electric field (fig. 7 and 8). Only by using the PULSTRODE PLUS catheters, the operator can also help himself in trying to hire the best treatment area, also using the dedicated switch integrated into the device.
- Only for PULSTRODE PLUS catheters, therefore select the size of the bipolar electric field most suitable for the treatment, choosing between the "big" position (proximal and distal active electrodes) and the "small" position (median and distal active electrodes). The "mono" position will instead exclude proximal and median electrodes and only work with the distal electrode (fig.7). Although less clinically effectiveness, with all catheters it is also possible to have a monopolar type of stimulation/treatment, which therefore provides the only active distal electrode and the application of a ground plate on the patient to be connected to the special plug socket on the external generator.


Fig. 7: electric field with 3 electrodes


Fig. 8: electric field with 2 electrodes

Pulstrode PLUS

EasyTrove PLUS

- At the end of the treatment remove the catheter and introducer, then medicate and put the patient at rest.

### RESPECT CAREFULLY

- Use the device only by qualified medical personnel.
- Observe instructions for use and those of the manufacturer of the Radio Frequency Generator.
- Ensure that the generator is in perfect working condition and has performed normal maintenance checks.
- Set the maximum working temperature limit of the generator @ 42° Celsius degree. Above this temperature the generator must automatically stop the signal delivery.
- Always check that the generator is set for PULSE RADIO FREQUENCY (PRF) treatments and NEVER for thermal radiofrequency treatments.
- NEVER perform treatments in temperature mode greater than 42° Celsius degree.
- NEVER use other introducers than the one included in the package.
- NEVER inject anesthetic in the epidural space before the treatment.
- NEVER run an angle bent on the catheter.
- NEVER perform treatments if the patient complains of unbearable burns or annoyances.
- NEVER infuse substances during treatment.
- NEVER perform treatments exceeding 75 Volt.
- NEVER modify the catheter position during the PRF procedure.
- NEVER place the catheter inside the subarachnoid space.
- NEVER turn off the generator with the catheter still inside the patient.

### COMPLICATIONS

It has been shown that catheters are safe and effective devices for the intended use. As in any surgical procedure of this type, if not properly performed, a number of complications could occur, such as increased pain, meningitis, epidural hemorrhage, dorsal puncture, continuous numbness, tingling, dysesthesia or paresthesia, arachnoiditis, subcutaneous infection or epidural.

### POSSIBLE COLLATERAL EFFECTS

In a percentage of patients after treatment, a recurrence of pain may occur for a period not exceeding two weeks.

### WARNINGS

In order to avoid:

- Infections: apply the best possible aseptic conditions and perform common perioperative antibiotic prophylaxis.
- Internal hemorrhages and hematomas during placement in a difficult area to reach: change the epidural access level and always maintaining a slight approach of the catheter when moving it through the epidural space.
- Thermal damage: set the generator with a safety max limit of temperature threshold and never exceeding 42° Celsius degree.
- Allergic reactions due to injected substances: observe the indications of the used medicaments.
- Structural catheter lesions: avoid insisting on penetrating too rigid scarring structures, adhesions or stenoses that obstruct the epidural canal. Avoid epidural access with the needle near the vertical of the patient's coronal plane (recommended angle is between 30° to 45° degree).

### REQUIRED EQUIPMENT FOR PULSE RADIOFREQUENCY TREATMENTS

For the correct PRF procedure, is required the following equipment:

- Pulsed radiofrequency generator with temperature monitoring and output signal parameters adjustable according to the following ranges:
  - PRF application time: adjustable with at least 300 sec reachable;
  - Treatment temperature: not higher than 42° Celsius degree;
  - PRF pulse width: adjustable and at least between 10 and 50 ms;
  - PRF pulse rate: adjustable and between 1 and 5 Hz;
  - Voltage treatment: adjustable with at least 75 volts reachable;
- Suitable generators: Cosman, RFG-1A model and newer;
- Connection cable for generator;
- Plate and connecting cable (only for monopolar use).

### AUTHORIZED USERS

The device can only be used by specialized medical staff and trained for this procedure, and only in the health sector both public and private. It is forbidden any other use in areas other than those mentioned above.

### RECOMMENDATIONS

Open the sterile bag immediately before the use to ensure its sterility. Do not use if the single package is open or damaged. Do not reuse the device: risk of onset of bacterial and viral infections.

### DESCRIPTION OF THE LABEL SYMBOLS

|  |  |  |   |
|--|--|--|---|
|  | REFERENCE CODE: This symbol is accompanied by the product identification code.   |  | MANUFACTURER: Indicates the name and address of the device manufacturer.  |
|  | LOT NUMBER: This symbol is accompanied by the manufacturer's batch code. Synonyms: "batch code" - "batch number".  |  | LATEX FREE: Indicates that the product and its primary packaging are latex-free.  |
|  | UNIQUE DEVICE IDENTIFIER: Indicates a vector containing information about the unique identification of the device.   |  | NO URBAN WASTE: Identifies that the product should not be disposed of as an urban waste, but separate collection should be made for it.       |
|  | MEDICAL DEVICE: Indicates that the item contained in the package is a medical device.  |  | DAMAGED PACKAGING: Do not use the product if the packaging is not intact.   |
|  | INSTRUCTIONS FOR USE: Indicates the need for the user to read the operating instructions before using the product for proper use.                            |  | KEEP AWAY FROM HUMIDITY: Indicates a medical device that must be protected from humidity.   |
|  | WARNING: It is placed as an alarm for special warnings to be observed before, during and after use of the product.   |  | HUMIDITY: Indicates the relative humidity range at which the product can be stored maintaining its features and functionality unchanged.      |
|  | STERILE PRODUCT: Sterilization at ETHYLENE OXIDE.  |  | STORAGE TEMPERATURE: Indicates the temperature range at which the product can be stored maintaining its features and functionality unchanged. |
|  | EXPIRY DATE: Date, expressed with four digits for the year and two digits for the month, indicating that the product should be used by the end of the month. |  | KEEP OUT OF SUNLIGHT: Indicates a medical device that must be protected from light sources.   |
|  | PRODUCTION DATE: Date, expressed with four digits for the year and two digits for the month, indicating the period of production.                            |  | DO NOT RESTERILIZE: Identify that the product can not be resterilized.  |
|  | DO NOT USE: This symbol identifies that the product is not reusable. CAUTION: Reuse can cause infections.  |  | IDENTIFIES THE NOTIFIER BODY: certifies the product's compliance with European Union regulations.   |

### PACKAGING

The device is sterilized with ethylene oxide gas and packed in a double tyvek pouch. The intact package guarantees sterility for 10 years. All catheters are disposable and not designed for resterilization and reuse.

### IN THE BOX FOR PULSTRODE PLUS (REF. PT46-4) AND EASYTRODE PLUS (REF. ET46-4)

The package includes:

- n° 1 electrocatheter with dispenser tube and stylet already inserted;
- n° 1 additional stylet complete with dispenser tube;
- n° 1 14-gauge steel introducer needle;
- n° 1 Torquer;
- n° 1 low resistance 10ml syringe;
- n° 1 cable cover tube;
- n° 1 blister tray with cover;
- n° 1 adhesive labels sheet with device info;
- n° 1 leaflet warnings for proper use.

### NOTES FOR STORAGE AND TRANSPORTATION

Environmental storage and transport conditions:

- temperature: 10° to 40° Celsius degree;
- relative humidity: 30% to 60%.

Store in dry areas, away from direct heat sources. Do not use the device after the expiration date. Keep out of the reach of children.

### DISPOSAL NOTES

If the device is not being used, then it must be disposed of separately waste as WEEE directive (Electrical and Electronic Equipment Rejection), vice versa if used, after use should be disposed as contaminated hospital waste. The penalties provided for the abusive disposal of such waste are set out in the regulations in actually force.

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BIOAMPERE RESEARCH Srl